

APPLICATION NOTE 2

¹Steinbeis Center for Biopolymer Analysis & Biomedical Mass Spectrometry,
Rüsselsheim am Main, Germany

and

²Reichert- Ametek Technologies, Buffalo, N.Y., USA

Epitope Identification and Affinity Determination of A β - specific Antibodies by online SPR- MS

Loredana Lupu¹, Hendrik Rusche¹, Zdenek Kukacka¹, Yannick Baschung¹, Mary Murphy²,
Jeff Bornheim², and Michael Przybylski¹

¹Steinbeis Center for Biopolymer Analysis & Biomedical Mass Spectrometry, Marktstraße 29,
Rüsselsheim am Main, Germany

²Reichert- Ametek Technologies, Walden Avenue, Buffalo, N.Y., USA

The accumulation of extracellular plaques containing the neurotoxic β -amyloid peptide fragment, A β (1-42) of β -amyloid precursor protein (β APP), is one of the characteristics of Alzheimer's disease (AD). Although β APP has been recognised as a key molecule for AD, its molecular (patho)physiological degradation, proteolytic pathways and cellular interactions of A β are still unclear. Studies towards the development of immunotherapeutic methods for AD have yielded initial success in transgenic mouse models of AD, in producing therapeutic antibodies by immunisation with A β (1-42), that disaggregate A β -plaques and fibrils. Using proteolytic excision of the immobilised A β antigen-immune complex in combination with ESI mass spectrometry, the A β -plaque specific epitope was identified as an N-terminal peptide A β (4-10), accessible in A β (1-42) as well as in oligomeric A β -fibrils [1, 2].

Recently, we have identified the epitope recognised by A β -autoantibodies in serum, capable of eliciting a neuroprotective effect to inhibit the formation of A β -plaques, to be located in the carboxyterminal region of the A β sequence. The differential epitope structures of A β -specific antibodies from healthy individuals and AD patients provides a breakthrough and molecular basis for (i), the development of new immunotherapeutic approaches by passive immunisation with A β -specific antibodies, and (ii), the development of new diagnostic tools for AD with absolute specificity (Figures 1; 2) [2, 3]. The primary structures of polyclonal A β -autoantibodies were elucidated by two-fold A β - epitope specific affinity chromatography from human immunoglobulin G, using a combination of overlapping proteolytic digestion (trypsin, α -chymotrypsin), HPLC isolation, and high resolution MS, which provided sequence data for Fv domains, CDR motifs, and framework regions. The plaque-protective A β -epitope of the A β -autoantibody was identified by proteolytic extraction-MS using the online SPRMS- epitope analyser, to reside in the A β (17-28) tryptic peptide sequence. To identify the epitope, the A β - autoantibody was immobilized on a dextran-SPR

affinity chip, and a tryptic mixture of A β -peptide fragments injected through the SPR autosampler followed by ESI-MS. Subsequent online desalting of analyte prior to MS was performed after elution of affinity captured A β -peptide which provided identification of the A β (17-28) epitope peptide (protonated molecular mass 1324,8). Following elution of undigested A β (1-40) through the microfluidic interface, SPR affinity determination revealed high affinity with a K_D of ca. 3.5 nM (Figures 3 and 4). Based on this epitope peptide, interactions of the plaque-protective A β -epitope with two fibril inhibiting peptides, cystatin-C and humanin, were evaluated at the molecular level to gain insight into the mode of action of A β -autoantibodies [4].

In this study we show that the online SPR- ESIMS combination is a powerful tool to enable the simultaneous affinity isolation, structure identification and affinity quantification of an A β - plaque protective epitope from the complex of A β -autoantibodies immobilized on a gold chip. The high application potential of online-SPR-MS has become further evident in recent studies of the identification of an unusual mixed-disulfide antibody epitope of the rheumatic target protein, HLA-B27; and the interaction site identification of chaperone complexes of lysosomal enzymes [5, 6]. Current applications confirm that interaction epitopes as diverse as antigen-antibody and lectin- carbohydrate complexes [7], and binding constants (K_D) from milli- to nanomolar ranges are amenable to SPR-MS analysis. These results indicate that applications of the online-SPR-MS epitope analyzer are well feasible to affinity-based biomarker evaluation; identification of protein and peptide epitopes; precise antibody affinity characterization; and direct label-free antigen quantification.

[1] McLaurin, J., R. Cecal, M.E. Kierstead, X. Tian, A.L. Phinney, M. Manea, J.E. French, M.H.L. Lambermon, A.A. Darabie, M.E. Brown, C. Janus, M.A. Chishti, P. Horne, D. Westaway, P.E. Fraser, H.T.J. Mount, M. Przybylski, P. St.-George-Hyslop (2002) *Nature Med.* **8**: 1263-1269.

[2] Stefanescu R., Iacob R.E., Damoc E.N., Marquardt A., Amstalden E., Manea M., Perdivara I., Maftai M., Paraschiv G., Przybylski M. (2007) *Eur. J. Mass Spectrom.* **13**: 69-75.

[3] Przybylski, M. et al./Univ. Konstanz & Bonn (2009) Eur. & US Patent Applications; M. Przybylski et al. (2011/2012), Eur. Patent Appl.

[4] Juszczak P, Paraschiv G, Szymanska A, Kolodziejczyk AS, Rodziewicz-Motowidlo S, Grzonka Z, Przybylski M. (2009) *J. Med. Chem.* **52**: 2420-2428.

[5] Iurascu, M.I., Marroquin Belaunzar, O., Cozma, C., Petrusch, U., Renner C., Przybylski M. (2016) *J. Am. Soc. Mass Spectrom.* **27**: 1105-1112.

[6] Moise, A., Maeser, S., Rawer, S., Eggers, F., Murphy, M., Bornheim, J., Przybylski, M. (2016) *J. Am. Soc. Mass Spectrom.* DOI: 1007/13361-016-1386-0.

[7] Moise A., Andre, S., Eggers, F., Krzeminski, M., Przybylski, M., Gabius, H.J. (2011) *J. Am. Chem. Soc.* **133**, 14844-14847.

Epitope identification and affinity determination of β -amyloid antibody




Fig 1: Abeta epitope peptides resulting from the cleavage of the molecule by APP. Abeta epitopes against antibodies are highlighted.




Fig. 2: Schematic representation of Amyloid protein binding to the A β - specific antibody. The SPR response increases exponentially during the binding and decreases upon epitope elution.

Epitope identification and affinity determination of β -amyloid antibody




Fig. 3: SPR determination of a dilutions series of $A\beta(1 - 40)$ upon processing via the SPR-MS interface. Kinetic evaluation resulted in a calculated K_D of 3.56 nM.




Fig. 4: ESI-MS identification by online SPR-MS of the epitope $A\beta(17-28)$ eluted from the $A\beta$ -antibody upon proteolytic extraction.